



Certificate

Herewith we, OMRON HEALTHCARE Co., Ltd. certify the blood pressure monitor HEM-6235 is substantially equivalent to HEM-6232T* which is clinically validated and employ the same algorithm as the predicate device for determining blood pressure.

Very Truly Yours,

*Validation of two automatic devices, Omron HEM-6232T and HEM-6181, for self-measurement of blood pressure at the wrist according to the ANSI/AAMI/ISO 81060-2:2013 protocol and the European Society of Hypertension International Protocol revision 2010
Kanao Saito , Yukiko Hishiki , Hakuo Takahashi

April, 28, 2021

Takehiro Hamaguchi
Technology Development HQ
Senior General Manager